

Tamper-evident: do not use if the imprinted seal under the cap is broken or missing.

PEEL

Expectorant Cough Syrup

Drug Facts

Active ingredient (in each teaspoonful (5 mL)) **Purpose**
Guaifenesin, USP 100 mg Expectorant

Uses

helps loosen phlegm (mucus) and thin bronchial secretions to:
 drain bronchial tubes
 make coughs more productive
 rid the bronchial passageways of bothersome mucus

Warnings

Ask a doctor before use if you have

cough that occurs with too much phlegm (mucus)
cough that lasts or is chronic such as occurs with smoking, asthma, chronic bronchitis, or emphysema ►

168x2.1

Outside label

Drug Facts (continued)

Stop use and ask a doctor if cough lasts for more than 7 days, comes back, or occurs with a fever, rash, or headache that lasts. These could be signs of a serious condition.

If pregnant or breast-feeding, ask a health professional before use.

Keep out of reach of children. In case of overdose, get medical help or contact a Poison Control Center right away.

Directions

take every 4 hours

do not take more than 6 doses in 24 hours

adults and children 12 years and over	2 to 4 teaspoonfuls (10 to 20 mL)
---------------------------------------	-----------------------------------

children 6 to under 12 years	1 to 2 teaspoonfuls (5 to 10 mL)
------------------------------	----------------------------------

children 2 to under 6 years	1/2 to 1 teaspoonful (2-1/2 to 5 mL)
-----------------------------	--------------------------------------

children under 2 years	ask a doctor ►
------------------------	----------------

Drug Facts (continued)

Other information

each teaspoonful (5 mL) contains: sodium 2 mg
store at room temperature 15 °C to 30 °C (59 °F to 86 °F)
protect from freezing protect from light
alcohol free
to report serious side effects due to the use of this product call (515) 276-1586

Inactive ingredients

caramel, citric acid, FD&C red no. 40, flavors, glucose, glycerin, high fructose corn syrup, menthol, sodium saccharin, sodium benzoate, water

Inside label

Base label - adheres to bottle